

**EXHIBIT I-2**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)**  
**ARIZONA DEPARTMENT OF PUBLIC SAFETY**  
**INTOXILYZER MODEL 9000**

**PERIODIC MAINTENANCE AND STANDARD QUALITY ASSURANCE PROCEDURE**

QA SPECIALIST \_\_\_\_\_ AGENCY \_\_\_\_\_  
DATE \_\_\_\_\_ TIME \_\_\_\_\_  
INTOXILYZER SERIAL # \_\_\_\_\_

1. Ensure that gas tank is attached and contains a standard alcohol concentration \_\_\_\_\_ AC

**DIAGNOSTIC TESTS**

1. Clock time check  
 2. Date check

**OPERATIONAL TESTS**

1. Deficient Subject Test (Proper Sample Recognition)  
Deficient Sample printed
2. Alcohol-free Subject Test (Proper Sample Recognition):  
0. \_\_\_\_\_ AC
3. Mouth Alcohol Subject Test (Proper Sample Recognition):  
Invalid Sample – Begin new deprivation period printed
4. Radio Frequency Interference Test (Error Recognition):  
RFI Detect printed
5. Standard Calibration Check:  
0. \_\_\_\_\_ AC
6. Air Blanks Completed
7. Timer Reset

Not a Successfully Completed Test Sequence will be printed

Instrument is operating properly and accurately. YES \_\_\_\_\_ NO \_\_\_\_\_

**COMMENTS:**

---

---

---

SIGNATURE \_\_\_\_\_