

**EXHIBIT G-2**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)**

**ARIZONA DEPARTMENT OF PUBLIC SAFETY**  
**STANDARD QUALITY ASSURANCE PROCEDURES**  
**INTOXILYZER MODEL 8000**

**STANDARD CALIBRATION CHECK PROCEDURE**

QA SPECIALIST \_\_\_\_\_ AGENCY \_\_\_\_\_

DATE \_\_\_\_\_ TIME \_\_\_\_\_

INTOXILYZER SERIAL # \_\_\_\_\_ LOCATION \_\_\_\_\_

- ( ) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution \_\_\_\_\_AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ( ) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- ( ) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- ( ) 4. Air blank completed.
- ( ) 5. Calibration check completed. Test results 0. \_\_\_\_\_AC.
- ( ) 6. Air blank completed.
- ( ) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_