Agilent TPS-compact

User Manual

Notices

Manual Part Number

87-901-010-01
Edition L.00, April 2022

Copyright

© Agilent Technologies, Inc. 2022

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Agilent Technologies Italia S.p.A.
Vacuum Products Division
Via F.lli Varian, 54
10040 Leini (TO)
ITALY

www.agilent.com

Instrument Manufacturing

Manufactured by Agilent Technologies Italia S.p.A., Vacuum Products Division
Via F.lli Varian, 54, 10040 Leini (TO),
ITALY

Printed in Italy

Operating Temperature

Operating Temperature: + 5 °C to + 35 °C
Storage Temperature: - 20 °C to + 60 °C

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data - Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

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Contents

7	Instructions for Use	157
	About this manual	158
	Validity	158
	Definitions and terminology	159
	Definition of Caution, Warning and Note	159
	Warning Symbols	160
	Safety	162
	Proper use	162
	Improper use	163
	Protective equipment	164
	Safety guideline for Turbomolecular Pumps	165
	Transport & Storage	166
	Product description	167
	Installation	168
	Preparation for installation	168
	Set-up	169
	TPS-compact Mounting	170
	Connection of the High Vacuum Flange	172
	Installation of ISO-K Flange	173
	Installation of ConFlat Flange	174
	Startup and Operation of the TPS-compact	175
	Soft Start	176
	Stopping the TPS-compact	177
	Emergency stop	177
	Maintenance	178
	Rotor life	178

Contents

Cleaning	178
Disposal	179
Service	180
8 Technical Information	181
Description of the TPS-compact	182
Turbopump Description	184
Controller Description	184
Control Panel Description	185
Scroll Pump Description	186
Technical specification	187
TPS-compact Outline	189
Inlet screen installation	192
High Vacuum Flange Connection Configurations	195
ASDF – Adaptive Scroll Driving Firmware	196
Gas Ballast	198
Interconnections	200
Input Power Connector	200
P3 – Vent	200
Gauge Connector	201
Pressure Display Module	202
J5 – In-Out	203
Signal Description	204
How to Connect the Open Collector Input of the Controller	209
How to Connect the Outputs of the Controller	212
J6 – Serial	212

Contents

RS 232/RS 485 Communication Description	214
Communication Format	214
Serial Communication Descriptions	216
Window Protocol	217
Description	217
Communication Format	217
Communication Protocol	217
Examples	220
Window Meanings	223
Service Options	226
Tip-Seal Replacement	227
To Remove the Worn Tip-Seals	228
To Install the new Seals and O-Rings	229
To test the Pumping System	230
Troubleshooting Chart	231
Accessories and Spare Parts	232

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7

Instructions for Use

About this manual	158
Validity	158
Definitions and terminology	159
Definition of Caution, Warning and Note	159
Warning Symbols	160
Safety	162
Proper use	162
Improper use	163
Protective equipment	164
Safety guideline for Turbomolecular Pumps	165
Transport & Storage	166
Product description	167
Installation	168
Preparation for installation	168
Set-up	169
TPS-compact Mounting	170
Connection of the High Vacuum Flange	172
Installation of ISO-K Flange	173
Installation of ConFlat Flange	174
Startup and Operation of the TPS-compact	175
Soft Start	176
Stopping the TPS-compact	177
Emergency stop	177
Maintenance	178
Rotor life	178
Cleaning	178
Disposal	179
Service	180

About this manual

Validity

This manual lists the instructions for the users of the TPS-compact, with particular reference to the notions relating to safety, operation and first level maintenance, limited to maintenance operations for which the user is responsible.

The maintenance operations, illustrated in the specific sections, with specific provisions relating to the higher level of maintenance (personnel specifically trained for maintenance operations) must not be carried out by the user.

For a correct installation and start/stop, please refer to "Installation" section.

Detailed information is supplied in the "Technical Information" section.

NOTE

- 1 This manual contains useful information so that all personnel using the TPS-compact can operate it safely and guarantee perfect efficiency, for its entire life span.
 - 2 Keep this manual, together with all the related publications, in an accessible place known to all operators/maintenance personnel.
-

Definitions and terminology

Definition of Caution, Warning and Note

Some important references of this manual are highlighted and framed in contrasting color.

CAUTION

Caution messages are displayed before procedures which, if not observed, could cause damage to the equipment.

WARNING

Warning messages draw the operator's attention to a specific procedure or practice which, if not performed correctly, could result in serious personal injury.

NOTE

Notes are intended to call attention to important information and provide more detail regarding specific steps.

Warning Symbols

The following is a list of symbols that appear in conjunction with warnings on the TPS-compact. The hazard they describe is also shown.

A triangular symbol indicates a warning. The meanings of the symbols that may appear alongside warnings in the documentation are as follows:

Protective Ground Terminal

Do Not Touch

European Declaration of Conformity

Generic hazard

Dangerous voltages

Hot surfaces

Manufacturing Site

Power OFF

Power ON

Restriction of Hazardous Substances certification

UK CA Marking

Instructions for Use

The following symbol may be used on warning labels attached to the instrument. When you see this symbol, refer to the relevant operation or service manual for the correct procedure referred to by that warning label.

The following symbols appear on the instrument for your information.

	Generic hazard
	CE certification
	CSA certification
	RoHS China certification
	Waste Electrical and Electronic Equipment
	UK CA Marking

Safety

This section contains the information, prescribed by the Machinery Directive 2006/42/EC, which is essential for the compliance and observance of the safety regulations both generally and in relation to the specific use of the machine.

Failure to comply with these instructions and the other instructions contained in this manual may render the safety conditions envisaged in the design phase inefficient and cause accidents to those operating the machine.

Agilent Technologies declines all responsibility for damage to the machine or for the physical safety of the operator or third parties deriving from the non-observance of the safety rules indicated in the technical documentation.

Proper use

This manual contains important warnings and safety instructions to be observed in order for the unit to work safely.

The product described in this manual is intended exclusively for the area of application specified in the instructions. The manual also provides indications regarding the essential requirements for the application and operation of the product as well as the safety measures that can be adopted to guarantee regular operation. Agilent Technologies does not provide any guarantee or assume any responsibility for applications other than those described in this manual or in which the essential requirements and safety measures are not respected.

The product must only be used by qualified personnel who are able to take the necessary safety measures under conditions that do not cause damage or injury. Any accessories and equipment used with the product must be supplied or approved by Agilent Technologies.

Any adjustment or maintenance operation must be performed by a professional technician informed about the risks.

Repairs on the product must be carried out exclusively by Agilent authorized personnel.

Improper use

Agilent Technologies declines all responsibility, deriving from the improper use of the TPS-compact.

Improper use will cause all claims for liability and warranties to be forfeited. Improper use is defined as:

- installation of the pump with unspecified mounting material
- pumping corrosive and oxidizing gases
- pumping of explosive gases
- pumping of condensing vapors
- pumping of liquids
- pumping dust
- operation with improper high gas throughput
- operation with improper high fore-vacuum pressures
- operation with selected gas modality different from the gas actually used
- operation with high levels of heat input
- operation in improper high magnetic fields
- venting with improper high venting rates
- use of the vacuum pump to generate pressure
- operation of the device in areas with ionizing radiation
- operation in potentially explosive areas
- using the system in an environment where stresses and vibrations lead to periodically reoccurring shocks and forces.

Protective equipment

The protective equipment of the operators who are operating or executing the maintenance of the pumping system must always be adequate for the type of operation being executed. Furthermore, it must comply with the safety requirements of the legislation in force in the country in which the machine is used.

In general, the operator must wear accident-prevention shoes while handling the TPS-compact and during installation.

WARNING

Danger to health by hazardous substances during maintenance or installation.

- Depending on the process peculiarity, vacuum pumps, components or operating fluids can be contaminated by toxic, reactive or radioactive substances.
- Wear adequate protective equipment during maintenance and repairs or in case of reinstallation.

WARNING

Risk of injury through falling objects

When transporting vacuum pumps by hand, there is a danger through loads slipping and falling down.

- Carry small and mid-size vacuum pumps two-handed.
- Any equipment heavier than 20 kg should be transported using a suitable lifting device.
- Wear safety shoes with steel toe cap according to directive EN 347.

Instructions for Use

CAUTION

Risk of injury through sharp edges

- Before servicing the pump or before any mounting/ dismantling action of the pump from the system, wait for the complete standstill of the pump.
- Do not operate directly inside the high vacuum flange.
- If necessary, wear protective gloves according to EN 420.

Safety guideline for Turbomolecular Pumps

Turbomolecular pumps as described in the following operating manual contain a large amount of kinetic energy due to the high Rotation Speed in combination with the specific mass of their rotors.

In case of malfunctioning of a mechanical part of the pump the rotational energy may be released to whatever the pump is fixed to (see "Installation of ConFlat Flange" and "Installation of ISO-K Flange" sections).

WARNING

To avoid damage to equipment and to prevent injuries to operating personnel the installation instructions as given in this manual should be strictly followed!

WARNING

When employing the pump for pumping toxic, flammable, or radioactive gases, please follow the required procedures for each gas disposal. The pump is designed to pump high throughput of N₂, Ar and lighter gas. Should you need to pump gases heavier than Ar, please contact Agilent technical support for informations.

The TPS-compact is not intended to be used in presence of aggressive, corrosive or explosive gases.

Transport & Storage

In order to guarantee the maximum level of performance and reliability of Agilent Turbomolecular pumps, the following guidelines must be followed to avoid condensation.

- when shipping, moving and storing pumps, the following environmental specifications should not be exceeded:
 - temperature range: refer to the manuals of the single TPS-compact component.
 - relative humidity range: refer to the manuals of the single TPS-compact component.
- the turbomolecular pumps must be always soft-started when received and operated for the first time by the customer.
- the turbomolecular pumps must be always soft-started after a period of 60 days of inactivity of the pump.

CAUTION

If for any reason the shelf life time is exceeded, the pump has to be returned to the factory. Please contact the local Agilent Vacuum Sales and Service representative for information.

WARNING

Use extreme caution handling heavy parts. A two person lift is recommended.

Failure to perform a two person lift may result in personal injury.

Product description

This equipment is intended for professional use and to be used by professionals.

The TPS-compact is an integrated system with a turbomolecular pump (TwisTorr 74 FS, Twistorr 305 FS) for high and ultra-high vacuum applications associated with its relevant controller and its primary pump. The system can pump any type of gas or gas compound. It is not suitable for pumping liquids or solid particles.

The pumping action is obtained through a high speed turbine (max. 70000 rpm) driven by a high-performance electric motor.

It is free of contaminating agents and, therefore, is suitable for applications requiring a "clean" vacuum.

The TPS-compact is equipped with auxiliary connectors to be controlled from a remote site by means of an host computer connected through a serial line.

Detailed information is supplied in the "Technical Information" section.

Installation

Preparation for installation

The TPS-compact is supplied in a special protective packaging. If this shows signs of damage which may have occurred during transport, contact your local sales office. Avoid sudden impacts to the TPS-compact that would be caused by dropping it. Be aware to maintain the packaging box in order to reuse it in case of storage for a long period of time or in case of shipping.

Do not dispose of the packing materials in an unauthorized manner. The material is 100% recyclable and complies with EEC Directive 94/62 and subsequent amendments.

CAUTION

In order to prevent outgassing problems, do not use bare hands to handle components which will be exposed to vacuum. Always use gloves or other appropriate protection.


Figure 31

NOTE

Normal exposure to the environment cannot damage the TPS-compact. Nevertheless, it is advisable to keep it closed until it is installed in the system, thus preventing any form of pollution by dust.

Set-up

CAUTION

Do not remove the adhesive and protective cap until you are ready to install the turbopump to the system.


Figure 32

Install this device indoors only, and do not use the system in an environment exposed to atmospheric agents (rain, ice, snow), dust, aggressive gases, or in explosive environments or those with a high fire risk.

During operation, the following environmental conditions must be complied with in order to avoid condensation:

- maximum ambient pressure: 2 bar above atmospheric pressure
- temperature: refer to the manuals of the single TPS-compact components
- relative humidity: refer to the manuals of the single TPS-compact components.

In the presence of magnetic fields the pumping system must be protected using a ferromagnetic shield.

TPS-compact Mounting

WARNING

Failure to comply with these installation instructions could result in the pump detaching from the system in the event of a rotor failure, which could cause property damage or serious injury or death.

The TPS-compact can be installed in vertical position only.

Fix the TPS-compact in a stable position, mounting the inlet flange of the turbopump to the system counter-flange.

The turbopump with ConFlat and ISO K inlet flange must be fixed to the vacuum chamber by means of the appropriate Agilent hardware.

For further details see the “Installation of ConFlat and ISO K” section.

Power supply cord: The correct cable for electrical wiring is a three wires (Ph+N+Earth) cable. The wire section has to be at least 0.75 mm² (AWG18).

CAUTION

The TPS-compact cannot be fixed by means of its base or through other flanges that are not inlet.

WARNING

The System must be installed in a way that allows an easy interruption of the line voltage.

CAUTION

To allow proper air circulation assure a free space all around the system at minimum of 20 cm and do not cover or reduce the space between the base of TPS-compact and its support plane.

Instructions for Use

For installation of optional accessories, see "Technical Information".

This paragraph details the fundamental operating procedures.

Make all electrical pneumatic connections before the use of the system. While heating the vacuum chamber, the temperature of the inlet flange must not exceed 120 °C with ConFlat flanging and 80 °C with ISO flanging.

WARNING

Never use the pumping system when the turbo inlet flange is not connected to the vacuum chamber or is not blanked.

WARNING

Do not touch the turbopump or any of its accessories during the heating process. The high temperatures may cause burns.

WARNING

Avoid impacts or harsh movements of the pump when in operation. The bearings may become damaged and damages to the persons or the things could be taken place.

CAUTION!

Use inert gas free from dust, particles or humidity (like Nitrogen) for venting the pump. The pressure at the vent port must be less than 2 bar (above atmospheric pressure).

WARNING!

Do not use the pumping system in presence of explosive, toxic, flammable or radioactive gases.

Connection of the High Vacuum Flange

If the rotor stops suddenly, the resulting torque must be absorbed by the hardware fastening the high vacuum flange.

Bolts or hardware must meet the specification below.

WARNING

It is mandatory to connect the pump in such a way to withstand the torque specified in the "TPS-compact Mounting" section. Specifically pay attention to clamp design, material of clamps and bolts and bolt fixing torque.

Use only high resistance class clamps (or bolts) ≥ 8.8 .

WARNING

If the pump is fastened on a vacuum chamber with flanges different from those recommended by Agilent, there is the chance of the TPS-compact twisting or tearing-off in case the rotor suddenly stops. The combination of different flange types is not allowed. Agilent will not accept any liability for all damages resulting from improper fastening.

Installation of ISO-K Flange

For ISO-K flange connections, fix the two flanges with the clamps model IC63250DCMZ or equivalent, as shown in the following figure (a protective screen can optionally be used).


Figure 33

Use the required number of claw clamps:

- 4 for ISO 63 flange
- 4 for ISO 100 flange
- 4 for ISO 160 flange

Tighten the bolts with a torque of 22 Nm.

Installation of ConFlat Flange

For ConFlat flange connections, fix the two flanges with the screws as shown in the following figure (a protective screen can optionally be used).


Figure 34

Use the required number of screws:

- 16 for CFF 6"
- 20 for CFF 8"

Attach the units and tighten each screw in turn. Repeat the sequential tightening until the flange faces meet. Tighten the screws with a torque of 20 Nm.

For ConFlat flange we recommend the usage of high strength steel bolts only.

To facilitate assembly and dismantling, it is optional to apply

Felpro C-100 high temperature lubricant to the screw threads protruding from the flange and between the nuts and flange.

In some cases, the connections can be made only with the bolt in the lower side.

CAUTION!

Be careful when tightening nuts and bolts to avoid creating dents on the turbomolecular pump casing.

Startup and Operation of the TPS-compact

Check the system voltage on the Agilent label stuck on the system electronic units case.

To switch on the TPS-compact it is sufficient to supply the mains and then move the external switch to ON position. The integrated controller automatically recognizes the mains presence and start up the pump.

WARNING

The initial startup of the pump TwisTorr 305 FS will be in "SOFT START" mode; for detailed information on this mode refer to the dedicated section.

The blue LED "STATUS" located on the TPS-compact base rear panel indicates with its flashing frequency the system operating conditions:

- with no flashing: the pump is normally rotating;
- slowly flashing (period of about 400 ms): the system is in ramp, or in braking, or in Stop;
- fast flashing (period of about 200 ms): error condition.

See the appendix "Technical Information" for a detailed description of the TPS-compact control panel.

Soft Start

WARNING

To ensure proper lubrication of the bearings the pump, at the first ramp-up, the pump has to be started in “soft start” mode.

This starting modality (Soft Start) ensures a smooth up to speed of the pump and it allows a proper lubrication of the bearings.

The control unit which is connected to the pump is already set-up by default to run the pump (only at the first control unit start-up). In case of use of a new pump with an used control unit, remember to set manually the soft start mode.

The soft start mode will take a longer time compared to the standard start-up; the user must wait that the pump is completely up to speed before stopping the pump.

After the first start, the pump won't be started in soft start mode again by the control unit.

NOTE

After 60 days or more of inactivity, in order to allow again the correct redistribution of the grease present in the bearings, it is necessary to start the pump via soft start mode again.

WARNING

Please refer to the controller's relevant instruction manuals in order to find the appropriate set of instructions to activate the “soft start” modality (Remote: Chapter Technical Information – Controller Connection & Window Meanings ; On Board: Chapter Technical Information – Interconnections & Window Meanings).

Stopping the TPS-compact

To switch off the TPS-compact it is sufficient to move the external switch to OFF position. The integrated controller immediately stops the pumping system.

WARNING

Before carrying out any work on the system, disconnect it from the power supply, vent the pump with the vent valve, wait until the rotor is completely stopped and wait until the surface temperature of the pump is below 50° C.

Emergency stop

To immediately stop the TPS-compact in an emergency condition it is necessary to remove the supply cable from the mains plug or moving the external switch to OFF position.

WARNING

This mode is to be used only in an emergency. The use of this stop mode could lead to faults and/or damages at the pump. Improper use of this stop mode could void the warranty.

Maintenance

The TPS-compact does not require any maintenance (except Tip-seal replacement). Any work must be carried out by authorized personnel.

In the case of breakdown, contact your local Agilent service center.

If a system needs to be scrapped, ensure that local regulations are followed for appropriate disposal.

Rotor life

The fatigue resistance of the turbomolecular pump rotor strongly depends on pump size. Refer to the specific pump manual for correct value (number of cycles of ramp up from stop to full speed and back to stop). As a precautionary measure, considering the typical usage of the pump and taking into account safety limits, Agilent recommends to return the system after 17 years of use for maintenance.

Cleaning

The exterior surface of the TPS-compact may be cleaned with mild detergent only.

Disposal

Meaning of the "WEEE" logo found in labels.

The following symbol is applied in accordance with the EC WEEE (Waste Electrical and Electronic Equipment) Directive.

This symbol (valid only in countries of the European Community) indicates that the product it applies to must NOT be disposed of together with ordinary domestic or industrial waste but must be sent to a differentiated waste collection system. The end user is therefore invited to contact the supplier of the device, whether the Parent Company or a retailer, to initiate the collection and disposal process after checking the contractual terms and conditions of sale.

Figure 35 Logo "WEEE"

For more information refer to:

<http://www.agilent.com/environment/product/index.shtml>

Service

Should a customer need an advanced exchange or repair service, please contact local distributor or directly mail to

vpt-customer@agilent.com

vpl-customer@agilent.com

Completion of the "Request for Return" form is required to return your pump to Agilent for service (provided at the end of this manual).

8 Technical Information

Description of the TPS-compact	182
Turbopump Description	184
Controller Description	184
Control Panel Description	185
Scroll Pump Description	186
Technical specification	187
TPS-compact Outline	189
Inlet screen installation	192
High Vacuum Flange Connection Configurations	195
ASDF – Adaptive Scroll Driving Firmware	196
Gas Ballast	198
Interconnections	200
Input Power Connector	200
P3 – Vent	200
Gauge Connector	201
Pressure Display Module	202
J5 – In-Out	203
Signal Description	204
How to Connect the Open Collector Input of the Controller	209
How to Connect the Outputs of the Controller	212
J6 – Serial	212
RS 232/RS 485 Communication Description	214
Communication Format	214
Serial Communication Descriptions	216
Window Protocol	217
Description	217
Communication Format	217
Communication Protocol	217
Examples	220
Window Meanings	223
Service Options	226
Tip-Seal Replacement	227
To Remove the Worn Tip-Seals	228
To Install the new Seals and O-Rings	229
To test the Pumping System	230
Troubleshooting Chart	231
Accessories and Spare Parts	232

Description of the TPS-compact

The TPS-compact pumping system consists of a pump with an integrated controller and a forepump and is available in twenty-four models which differ in the TPM size, high vacuum flange and voltage. The models are:

Table 1

P/N	TMP	Flange	Foreline	Voltage
X3583-64000	74 FS	KF40	3.6 m3/h Scroll	115V 60Hz
X3583-64001	74 FS	CFF 2.75	3.6 m3/h Scroll	115V 60Hz
X3583-64002	74 FS	ISO63	3.6 m3/h Scroll	115V 60Hz
X3583-64003	74 FS	CFF 4.5	3.6 m3/h Scroll	115V 60Hz
X3583-64004	74 FS	KF40	3.0m3/h Scroll	230V 50-60Hz
X3583-64005	74 FS	CFF 2.75	3.0m3/h Scroll	230V 50-60Hz
X3583-64006	74 FS	ISO63	3.0m3/h Scroll	230V 50-60Hz
X3583-64007	74 FS	CFF 4.5	3.0m3/h Scroll	230V 50-60Hz
X3583-64008	74 FS	KF40	3.6 m3/h Scroll	100V 50-60Hz
X3583-64009	74 FS	CFF 2.75	3.6 m3/h Scroll	100V 50-60Hz
X3583-64010	74 FS	ISO63	3.6 m3/h Scroll	100V 50-60Hz
X3583-64011	74 FS	CFF 4.5	3.6 m3/h Scroll	100V 50-60Hz
X3582-64000	305 FS	ISO100	3.6 m3/h Scroll	115V 60Hz
X3582-64001	305 FS	CFF6	3.6 m3/h Scroll	115V 60Hz
X3582-64002	305 FS	ISO160	3.6 m3/h Scroll	115V 60Hz
X3582-64003	305 FS	CFF8	3.6 m3/h Scroll	115V 60Hz
X3582-64004	305 FS	ISO100	3.0m3/h Scroll	230V 50-60Hz
X3582-64005	305 FS	CFF6	3.0m3/h Scroll	230V 50-60Hz
X3582-64006	305 FS	ISO160	3.0m3/h Scroll	230V 50-60Hz
X3582-64007	305 FS	CFF8	3.0m3/h Scroll	230V 50-60Hz
X3582-64008	305 FS	ISO100	3.6 m3/h Scroll	100V 50-60Hz
X3582-64009	305 FS	CFF6	3.6 m3/h Scroll	100V 50-60Hz
X3582-64010	305 FS	ISO160	3.6 m3/h Scroll	100V 50-60Hz
X3582-64011	305 FS	CFF8	3.6 m3/h Scroll	100V 50-60Hz

Technical Information


Figure 36


Figure 37

Technical Information

Turbopump Description

The turbopump consists of a high frequency motor driving a turbine fitted with many bladed stages and Twistorr drag stages. The turbine rotates in an anticlockwise direction when viewed from the high vacuum flange end.

The turbine is made of high-strength, light aluminium alloy, and is machined from a single block of aluminium. The Twistorr stages are in the form of discs.

The turbine rotor is supported by permanently lubricated high precision ceramic ball bearings installed on the forevacuum side of the pump.

The static blades of the stator are fabricated in stainless steel. These are supported and accurately positioned by spacer rings.

The TwisTorr stators are in the form of selfpositioning machined disks. On each disk there are parallel spiral pumping channels designed to pump in centrifugal direction on one side and in centripetal direction on the other side. These are made of aluminium alloy.

During normal operation, the motor functions with a power feed at 54 Vac three-phase at 1167 Hz (74 FS) or 1010 Hz (305 FS). To reduce losses during start-up to a minimum, the frequency increases according to a ramp with a higher initial voltage/frequency ratio.

A thermistor sensor is mounted near the upper bearing to prevent the pump from overheating. The pump is balanced after assembly with a low residual unbalance amplitude.

Controller Description

The integrated controller is a solid-state frequency converter which is driven by a single chip micro-computer and is composed of a PCB which includes a power supply with a 3-phase AC output, analog and input/output section, microprocessor and digital section. The controller recognizes the mains presence and converts the single phase AC mains supply into a 3-phase, medium frequency output which is required to power the pump.

The controller can be operated by a remote host computer via the serial connection.

Control Panel Description

The following picture shows the TPS-compact control panel; the associated table lists the available commands. For a detailed description see the following paragraphs.


Figure 38

Table 2

Control name	Description
J5 IN-OUT	Input-output connector to remote control the TPS-compact
J6 SERIAL	Serial input-output connector to control the TPS-compact via an RS 232 or RS 485 connection
STATUS	Blue LED to show the system operating condition
P3 VENT	Connector to control the optional vent valve
GAUGE	Connector to control external active gauge
ON-OFF switch (not shown in the picture)	It is located on the rear side of the system. It is the power supply switch.

Technical Information

Scroll Pump Description

The TPS-compact internal scroll creates vacuum using a simple dual scroll mechanism in which one of the nested scrolls orbits about the other, creating moving zones of captured gas.

Gas enters the scroll set at the perimeter and is displaced and compressed toward the center hub where it is exhausted.

Figure 39

CAUTION!

Scroll Pump exhaust connection: To avoid overheating the Scroll pump, do not restrict the exhaust flow with long lengths of small diameter tubing.

Technical specification

The following tables list the main technical data of the TPS-compact.

Table 3 TPS-compact Technical Specification

Characteristic	305 FS	74 FS
Pumping speed	N ₂ : 250 l/s He 255 l/s Ar 250 l/s	N ₂ : 48 l/s (with inlet screen) He: 44 l/s (with inlet screen) H ₂ : 32 l/s (with inlet screen)
Base pressure*	<1.0 ⁻⁹ mbar <7.5 ⁻¹⁰ Torr	1x10 ⁻⁹ mbar 7.5x10 ⁻¹⁰ Torr
Inlet flange	ISO 100, ISO 160, CFF 6, CFF 8	KF 40, ISO 63, CFF 2.75, CFF 4.5
Rotational speed	60600 rpm (1010 Hz driving frequency)	70000 rpm
Start-up time	350 sec (it can be longer in case of soft start active)	144 s
Operating position	Vertical	
Operating ambient temperature	+ 5 °C to + 35 °C	
Bakeout temperature	100 °C at inlet flange max (CF flange) 75 °C at inlet flange max (ISO flange)	120 °C at inlet flange max. (CF flange) 80 °C at inlet flange max. (ISO flange)
Noise level	55.4 dB(A) at 1 meter	51.5 dB (A) at 1 meter
Lubricant	permanent lubrication	
Power supply (Line type): Input voltage and frequency:	100 Vac 50/60 Hz 115 Vac 60 Hz 200–230 Vac 50 – 60 Hz	
Max input power:	650 VA	
Protection fuse Type Time Delayed	115 Vac: 5 A 230 Vac: 3 A 100 Vac: 5 A	115 Vac: 4 A 230 Vac: 2.5 A 100 Vac: 4 A
Compliance with:	EN 55011 (Class A Group 1) EN 61326-1 EN 61010-1 EN 1012-2 EN 12100	

* According to standard DIN 28 428, the base pressure is that measured in a leakfree test dome, 48 hours after the completion of test dome bake-out, with a Turbopump fitted with a ConFlat flange.

Technical Information

Table 3 TPS-compact Technical Specification - Continued


Characteristic	305 FS	74 FS
Indoor Use		Only
Max Altitude		2000m
Installation category		II
Pollution degree		2
Power cable		115 Vac: US plug 3 meters long 230 Vac: EU or UK plug 3 meters long
Serial communication		RS232 and RS485 interface
Storage temperature		-20°C to + 60°C
Weight kg (lbs)	20.1 (44.3)	16.7 (36.8)

NOTE


When the TPS-compact has been stored at a temperature less than 5 °C, wait until the TPS-compact has reached the above mentioned temperature before switch the system on.

TPS-compact Outline

The following figure shows the TPS-compact outline. Dimensions are in mm.

74 FS

305 FS

Dimensions: millimeters (inches)

Figure 40

Technical Information

TPS-compact equipped with 305 FS: Pumping speed for Nitrogen.


Figure 41 Graph of pumping speed vs inlet pressure at full speed (1010Hz driving freq.)

Technical Information


Figure 42 Graph of pumping speed vs inlet pressure

TPS-compact 305 FS: 15l pumpdown


Figure 43 Graph of pressure vs time (w/o soft start)*

* Measurement obtained using both capacitive and UHV gauges Shape of curve in the range 1.0E-3 mbar – 1.0E-5 mbar is strongly affected by the type of gauge used and its accuracy.

Technical Information


Figure 44 Graph of pressure vs time (with a pumpdown volume of 1 litre)

Inlet screen installation

Installing the inlet screen will reduce the pumping speed of heavy gases (i.e. Argon, Nitrogen) of about 25% while for lighter gases (i.e. Helium and Hydrogen) the amount is of about 10%.


Figure 45

Technical Information


Figure 46

The inlet screens mod. 9699300, 9699309, 9699302 and 9699303 prevent the blades of the pump from being damaged by debris greater than 0.7 mm diameter.

The inlet screen is fitted in the upper part of the pump, as shown in the figure.


Figure 47

Technical Information

The screen can be mounted on each pump.

The screen can be removed as shown in the following figure.


Figure 48

High Vacuum Flange Connection Configurations


Figure 49


Figure 50

ASDF – Adaptive Scroll Driving Firmware

TPS-compact implements a function that safeguards the life of the scroll pump Tip-seal. The function is described in the diagram below where P is the turbomolecular pump power.

Inlet gas load is monitored measuring the turbomolecular absorbed power. A power threshold P_{TMP}^{OFF} is defined and if the turbomolecular power is below the P_{TMP}^{OFF} , then the foreline separation electromagnetic valve is closed and integrated scroll pump is switched off. Since the turbomolecular pump temporarily works without a backing pump, the power necessary to maintain a stable inlet pressure tends to increase. When the absorbed power exceeds the threshold power P_{TMP}^{ON} , then the integrated scroll pump is switched on and the separation electromagnetic valve is re-opened. So, the ASDF is able to extend the scroll Tip-seal life.


NOTE

A slight variation in turbomolecular pump inlet pressure can be detected; this variation is a normal response of the TPS-compact to primary pump shutdown.

By disabling this function, the scroll remains in operation as long as the customer does not switch off the entire system.

By disabling this function, warranty does not cover the short tip seal life.

To disable this function, put a jumper between pins 3-9 and 4-15 on I/O connector of Turbo controller:

This signal is no more used for interlock, it's now dedicated to disable/enable the adaptive driving firmware.

Technical Information


Figure 51

Gas Ballast

The TPS-compact packaging includes a standard gas ballast kit (P/N IDP3GBKIT) that can be used to prevent water and other condensates from accumulating within the foreline pump. The standard configuration has a seal screw installed in the 1/4-20 thread in the top of the exhaust manifold. This configuration can be used for dry applications. When appreciable amounts of water, solvents, etc., are being pumped, remove the seal screw and replace it with the provided adapter and 1/8" NPT sintered filter plug.

WARNING

Before install gas ballast adapter remove the power supply cable.

- 1 Disconnect the TPS-compact from electrical power.
- 2 Unlock the scroll pump support by unscrewing the screws.


Figure 52

- 3 Lift the TPS-compact in order to access the Scroll pump site.


Figure 53

Technical Information

- 4 Unscrew seal screw and install provided gas ballast adapter and 1/8" NPT sintered filter plug.


Figure 54

- 5 Lock the scroll pump support by screwing the screws.

CAUTION

Be careful during pump support closing to prevent damage to system components (vacuum pipes and connectors).

NOTES

- During the operation of the gas ballast feature the noise may increase.
- There is not a predefined amount of time to let the gas ballast open: it depends on the amount of moisture trapped in the system.

CAUTION

Do not exceed 2.8 Nm of torque when installing the adapter into the manifold.

For more details about gas ballast refer to the manuals of the single TPS-compact components.

CAUTION

Use the provided adapter and 1/8" NPT sintered filter plug in order to avoid dust entrance.

CAUTION

The gas ballast must be sealed whenever pumping any gas not intended to be vented to atmosphere.

Interconnections

Input Power Connector

The following figure shows the input power connector.


Figure 55

P3 – Vent


Figure 56

This is a dedicated 24 Vdc connector to control the optional vent valve (available on request).

The vent valve can be driven by controller automatically or by serial line.

Technical Information


Figure 57

Gauge Connector

TPS-compact is equipped with new generation gauge reading card able to drive/read a Agilent Full Range Gauge FRG-700.

The FRG-700 Full Range Gauge is a combined technologies gauge (Inverted Magnetron plus Pirani Gauge).

FRG-700 measure from 5×10^{-9} mbar to atmosphere (3.8×10^{-9} Torr to atmosphere) and its temperature range goes from 5 °C to 55 °C [41 °-131 °F].

Pressure data is available on three independent reading channels:

- 6** On-board display (data expressed in mbar, Torr or Pa) – (see Pressure Display Module picture).
- 7** Serial communication line (Win 224).
- 8** Analogical output.

Technical Information

Pressure gauge can be connected/disconnected from the TPS-compact during normal operation


Figure 58 Gauge Connector

Pressure Display Module

Figure 59 Pressure display module


Figure 60 FRG-700 gauge dimension

Technical Information

J5 – In-Out


Figure 61

This connector carries all the input and output signals to remote control the TPS-compact.

It is a 15-pins D type connector; the available signals are detailed in the table, the following paragraphs describe the signal characteristics and use.

Table 4

Pin N.	Signal name	Input/Output
1	START/STOP (+)	IN
2	START/STOP (-)	IN
3	ADAPTIVE SCROLL DRIVING DISABLE (+)	IN
4	ADAPTIVE SCROLL DRIVING DISABLE (-)	IN
5	SPEED SETTING (+)	IN
6	SPEED SETTING (-)	IN
7	SOFT START (+)	IN
8	SOFT START (-)	IN
9	+24 Vdc I _{max} = 30 mA	OUT
10	NORMAL OUTPUT (relay contact) V _{max} = 125V I _{max} = 200 mA	OUT
11	PROGRAMMABLE SET POINT	OUT
12	NORMAL OUTPUT (relay contact)	OUT
13	FAULT OUTPUT	OUT
14	PROGRAMMABLE ANALOG SIGNAL (+)	OUT
15	<ul style="list-style-type: none">GROUNDPROGRAMMABLE ANALOG SIGNAL (-)	OUT

Signal Description

START/STOP: input signal to start or stop the pump. Without the supplied cover connector the turbomolecular pump automatically starts, if you want to stop the turbo-pump from outside, you have to connect the START/STOP (+) signal to the +24 Vdc pin and the START/STOP (-) signal to the GROUND pin (this is the default configuration of the cover connector).

ADAPTIVE SCROLL DRIVING DISABLE: Input signal to disable Adaptive Scroll Driving function. With the supplied mating connector the Adaptive Scroll Driving Disable (+) input is connected to the +24Vdc pin and the Adaptive Scroll Driving Disable (-) input is connected to GROUND pin; in this condition the Adaptive Scroll Driving function is disabled.

SOFT START: this input is used to provide a "soft start" to the pump; in this condition the ramp-up time could be up to 20 min.

SPEED SETTING: PWM input signal to set the pump speed. The PWM signal characteristics must be the following:

- frequency: 100 Hz +/-20 %
- amplitude: 24 V max
- duty cycle range: from 25 % to 75 %
- corresponding to a pump speed from 1100Hz to 1167Hz (see the following diagram).

Technical Information


Figure 62

PROGRAMMABLE ANALOG SIGNAL: this output signal is a voltage (from 0 to 10 Vdc) proportional to a reference quantity (frequency or power or Temperature or Pressure) set by the user. The de-fault setting is the frequency (see the following example diagram).


Figure 63

Technical Information


Figure 64


Figure 65


Figure 66

FAULT: this open collector output signal is ON when a system fault condition is detected.

PROGRAMMABLE SET POINT: this open collector output signal is enabled when the reference quantity chosen (frequency, power or time or pressure or normal) is higher than the set threshold. The signal can be "high level active" (that is the output is normally OFF and becomes ON when activated), or "low level active" (that is the output is normally ON and becomes OFF when activated). Moreover, if the reference quantity is the frequency or the power or pressure drawn, it is possible to set the hysteresis (in % of the threshold value) to avoid bouncing.

For example:

- reference quantity: frequency
- threshold: 867 Hz
- hysteresis: 2 %
- activation type: "high level"

Technical Information


Figure 67

It is possible to delay the set point checking for a programmable delay time.

The PROGRAMMABLE SET POINT signal has the following default settings:

- reference quantity: frequency
- threshold: 867 Hz
- hysteresis: 2 %
- activation type: high level
- delay time: 0 second

NOTE

The A-Plus Software (optional) allows the operator to set all the programmable feature.

How to Connect the Open Collector Input of the Controller

Here below there are the typical connections of the open collector input of TPS-compact to an external system. Two cases are considered:

- 1 the customer supplies the 24 Vdc.
- 2 the customer does not supply the 24 Vdc.

Please note that on the connector a 24 Vdc, 30 mA, a GROUND signal and the open collector pin are available.


Figure 68 Case 1

Technical Information


Figure 69 Case 2 with relay utilisation

Technical Information


Figure 70 Case 2 with transistor utilisation

Technical Information

How to Connect the Outputs of the Controller

The following figure shows a typical logic output connection (relay coil) but any other device may be connected e.g. a LED, a computer, etc., and the related simplified circuit of the controller. The figure example refers to the programmable set point signal on pins 11 and 9.


Figure 71 Typical output connection 1

J6 – Serial


Figure 72

This is a 9 pin D-type serial input/output connector to control via an RS 232 or RS 485 connection the TPS-compact.

Technical Information

Table 5

P/N	Signal name
1	+5 V (OUT)
2	TX (RS232)
3	RX (RS232)
4	Not used
5	GND
6	A + (RS485)
7	Not used
8	B - (RS485)
9	Not used

A serial communication kit with a serial cable and the A-Plus Software is available (optional).

RS 232/RS 485 Communication Description

Both the RS 232 and the RS 485 interfaces are available on the connector J6.

The communication protocol is the same (see the structure below), but only the RS 485 manages the address field. Therefore to enable the RS 485 is necessary to select the type of communication as well as the device address by means of the software.

Communication Format

- 8 data bit
- no parity
- 1 stop bit
- baud rate: 600/1200/2400/4800/9600 programmable.


Figure 73 RS 232 Connection

Technical Information


Figure 74 RS 485 Connection

Technical Information

The communication port mating connector is supplied with the RS232 PCB (AMP/Cannon or equivalent 15-pin "D" type male connector).

For example, the Transmit data signal from controller (pin 2) must be connected to the host computer's receive data line (pin 2) and vice versa. Consult the host computer's instruction manual for its serial port connections

NOTE

Agilent cannot guarantee compliance with FCC regulations for radiated emissions unless all external wiring is shielded, with the shield being terminated to the metal shroud on the 0-subconnector. The cable should be secured to the connector with screws.

Serial Communication Descriptions

This unit can communicate by two different protocols:

- protocol (old system) "letter"
- "Window" protocol (new system).

These two protocols can be used as well with 232 or 485 media.

NOTE

Please use "Window" protocol for new development.

Window Protocol

Description

Communication Format

- 8 data bit
- no parity
- 1 stop bit
- baud rate: 600/1200/2400/4800/9600 programmable

Communication Protocol

The communication protocol is a MASTER/SLAVE type where:

- Host = MASTER
- Controller = SLAVE

The communication is performed in the following way:

- 1** the host (MASTER) send a MESSAGE + CRC to the controller (SLAVE);
- 2** the controller answer with an ANSWER + CRC to the host.

The MESSAGE is a string with the following format:

<STX>+<ADDR>+<WIN>+<COM>+<DATA>+<ETX>+<CRC>

Where:

NOTE

When a data is indicated between two quotes ('...') it means that the indicated data is the corresponding ASCII character.

Technical Information

- <STX> (Start of transmission) = 0x02
- <ADDR> (Unit address) = 0x80 (for RS 232)
<ADDR> (Unit address) = 0x80 + device number (0 to 31)
(for RS 485)
- <WIN> (Window) = a string of 3 numeric character indicating the window number (from '000' to '999'); for the meaning of each window see the relevant paragraph.
- <COM> (Command) = 0x30 to read the window, 0x31 to write into the window
- <DATA> = an alphanumeric ASCII string with the data to be written into the window. In case of a reading command this field is not present.

The field length is variable according to the data type as per the following table:

Table 6

Data Type	Field Length	Valid Characters
Logic (L)	1	'0' = OFF '1' = ON
Numeric (N)	6	'0', '1', '2' ... '9' right justified with '0'
Alphanumeric (A)	10	from blank to '_' (ASCII)

- <ETX> (End of transmission) = 0x03
- <CRC> = XOR of all characters subsequent to <STX> and including the <ETX> terminator. The value is hexadecimal coded and indicated by two ASCII character.

The addressed SLAVE will respond with an ANSWER whose structure depends from the MESSAGE type. When the MESSAGE is a reading command, the SLAVE will respond transmitting a string with the same structure of the MESSAGE.

NOTE

Using the RS 485 interface, the message structure remains identical to the one used for the RS 232 interface, the only difference being that the value assigned to the ADDRESS <ADDR>.

Technical Information

The controller can answers with the following response types:

Table 7

Response Type	Response Length	Response Value	Description
Logic	1 byte	-	After a read instruction of a logic window
Numeric	6 bytes	-	After a read instruction of a numeric window
Alphanumeric	10 bytes	-	After a read instruction of an alphanumeric window
ACK	1 byte	(0x6)	The command execution has been successfully completed
NACK	1 byte	(0x15)	The command execution has been failed
Unknown Window	1 byte	(0x32)	The specified window in the command is not a valid window
Data Type Error	1 byte	(0x33)	The data type specified in the command (Logic, Numeric or Alphanumeric) is not accorded with the specified Window
Out of Range	1 byte	(0x34)	The value expressed during a write command is out of range value of the specified window
Win Disabled	1 byte	(0x35)	The specified window is Read Only or temporarily disabled (for example) you can't write the Soft Start when the Pump is running)

Examples

Command: START

Source: PC

Destination: Controller

02	80	30	30	30	31	31	03	42	33
STX	ADDR	WINDOW			WR	ON	ETX	CRC	

Source: Controller

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

Command: STOP

Source: PC

Destination: Controller

02	80	30	30	30	31	30	03	42	32
STX	ADDR	WINDOW			WR	OFF	ETX	CRC	

Source: Controller

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

Technical Information

Command: SOFT-START (ON)

Source: PC

Destination: Controller

02	80	31	30	30	31	31	03	42	32
STX	ADDR	WINDOW			WR	ON	ETX	CRC	

Source: Controller

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

Command: SOFT-START (OFF)

Source: PC

Destination: Controller

02	80	31	30	30	31	30	03	42	33
STX	ADDR	WINDOW			WR	OFF	ETX	CRC	

Source: Controller

Destination: PC

02	80	06	03	38	35
STX	ADDR	ACK	ETX	CRC	

Technical Information

Command: READ PUMP STATUS

Source: PC

Destination: Controller (with address = 3)

02	83	32	30	3	30	03	38	37
STX	ADDR	WINDOW			RD	ETX	CRC	

Source: Controller (with address = 3 in stop status)

Destination: PC

02	83	32	30	35	30	30	30	30	30	30	03	38	37
STX	ADDR	WINDOW			DATA (STATUS)						ETX	CRC	

Command: REAL SERIAL TYPE

Source: PC

Destination: Controller (with address = 3 in 485 mode)

02	83	35	30	34	30	03	38	31
STX	ADDR	WINDOW			RD	ETX	CRC	

Source: Controller

Destination: PC

02	83	35	30	34	30	31	03	42	30
STX	ADDR	WINDOW			RD	DATA	ETX	CRC	

Window Meanings

N.	Read/Write	Data Type	Description	Admitted Values
000	R/W	L	Start/Stop (in remote/ mode the window is a read only)	Start = 1 Stop = 0
001	R/W	L	Low Speed Activation (in remote/ mode the window is a read only)	No = 0 Yes = 1 (default = 0)
008	R/W	L	Remote or Serial configuration	Remote = 1 Serial = 0 (default = 1)
100	R/W	L	Soft Start (write only in Stop condition)	YES = 1 NO = 0
101	R/W	N	R1 Set Point type	0 = Frequency 1 = Power 2 = Time 3 = Normal 4 = Pressure (default = 0)
102	R/W	N	R1 Set Point valve (expressed in Hz, W or s)	(default = 867)
103	R/W	N	Set Point delay: time between the pump start and the set point check (seconds)	0 to 999999 (default = 0)
104	R/W	L	Set Point signal activation type: the signal can be "high level active" or "low level active"	0 = high level active 1 = low level active (default = 0)
105	R/W	N	Set point hysteresis (in % of value)	0 to 100 (default = 2)
107	R/W	L	Active Stop (write only in stop)	0 = NO 1 = YES (default = 0)
108	R/W	N	Baud rate	600 = 0 1200 = 1 2400 = 2 4800 = 3 9600 = 4 (default = 4)
111	R/W	L	Analog output type: output voltage signal proportional to frequency or power	0 = frequency 1 = power (default = 1) 2 = Pump temperature 3 = Pressure
117	R/W	N	Low Speed frequency (Hz)	1100 to "Maximum rotational frequency" (win 120) (default = 1100)

Technical Information

N.	Read/ Write	Data Type	Description	Admitted Values	
				For 74 FS	For 305 FS
120	R/W	N	Rotational frequency setting (Hz)	1100 to 1167 (default = 1167)	700 to 1010 (default = 1010)
122	R/W	L	Set vent valve on/off	0 = Closed 1 = Open (default = 1)	
125	R/W	L	Set the vent valve operation	Automatic = 0 (see note 1.) On command = 1 (see note 2.)	
126	R/W	N	Vent valve opening delay (expressed in 0.2 sec)	0 to 65535 (corresponding to 0 to 13107 sec) (default = 15)	
147	R/W	N	Vent open time	0 = infinite 1 bit = 0.2 sec 1 to 65535	
155	R	N	Power limit applied	155	
157	R/W	N	Gas load type Select the gas load to the pump	Ø = N ² 1 = Ar	
162	R/W	A	R1 Set Point Pressure Valve Valid if Win. 101 = 4 Format X.X E. XX Where X = 0 to 9 s = + or -		
163	R/W	N	Pressure unit of measure	0 = mBar 1 = Pa 2 = Torr	
167	R/W	A	Stop speed reading Activates / deactivates the pump speed reading after Stop command	0 = disable 1 = enable	
200	R	N	Pump current in mA dc		
201	R	N	Pump voltage in Vdc		
206	R	N	Error code	Bit description: see the following figure	
216	R	N	Controller Air Temperature (°C).		
224	R	A	Pressure reading Format = X.X E XX		
226	R	N	Rotation Frequency (rpm)		
300	R	N	Cycle time in minutes (zeroed by the reset command)	0 to 999999	
301	R	N	Cycle number (zeroed by the reset command)	0 to 9999	
302	R	N	Pump life in hours (zeroed by the reset command)	0 to 999999	
310 to 347			Reserved to Agilent service		

Technical Information

N.	Read/Write	Data Type	Description	Admitted Values
358	R/W	N	Tip-seal life (R) Tip-seal life reset (W)	0 = reset
406	R	A	Program Listing Code	"XXXXXXXXXX"
407	R	A	Parameter Listing Code	"XXXXXXXXXX"
500	Reserved to Agilent service			
503	R/W	N	RS 485 address	0 to 31 (default = 0)
504	R/W	L	Serial type select	0 = RS 232
724	R/W	N	Run Up time (minutes)	0.30 15
725	R/W	L	Run Up time control	On = 1 1 Off = 0

NOTE

- 1 Automatic means that when the controller stops, the vent valve is opened with a delay defined by window n. 126; when the controller starts, the vent valve is immediately closed.
- 2 On command means that the vent valve is opened or closed by means of window n. 122.
- 3 These commands remove the pump motor supply for a short time. If they are used continuously, they can cause the pump stopping.


Figure 75 Windows N. 206 Bit Description

Service Options

Agilent pumping systems will provide many years of trouble-free operation if the maintenance procedures are observed.

Pumping System cleaning and internal scroll pump tip seal replacement are recommended when Pumping System TMP adsorbed current (with blanked inlet flange) has risen to an unacceptably high level (>500 mA) or when the message "TIP-SEAL SERV" appears on the TPS-compact front side display.

The warning message appears after 5000 h of operation.

If your TPS-compact exhibits humming or grinding noises from the scroll pump, a major overhaul should be done by Agilent Vacuum Technologies.

The parts needed for Tip-Seal replacement on the internal scroll pump are available in the kit described in Orderable parts Table. This kit contains seals and O-rings, and can be obtained from your Agilent dealer.

Tip-Seal Replacement

WARNING

Before replace Tip-Seals remove the power supply cable.

The parts and tools required to replace Tip-Seals are listed below:

- Tip-Seal Replacement Kit P/N IDP3TS
- 4 mm Allen wrench
- Razor blade or side-cutting pliers
- Compressed air (optional)


Figure 76 Tip-Seal replacement - Figure A


Figure 77 Tip-Seal replacement - Figure B


Figure 78 Tip-Seal replacement - Figure C

To Remove the Worn Tip-Seals

Disconnect the TPS-compact from electrical power.

Unlock the scroll pump support by unscrewing the screws (A) (see the “Tip-Seal replacement - Figure A”)

- 1** Lift the TPS-compact in order to access the Scroll pump site (you’ll hear a “click” noise as soon as the upper part of TPS-compact is locked).
- 2** Remove (4) M5 socket head bolts (B).
- 3** Remove the outboard housing axially off the frame (C).
- 4** Remove and discard the worn tip seals (D) and the main O-ring.
- 5** If compressed air is available, blow any remaining seal debris off the scroll parts. If seal debris is attached to the sides, use a razor or Exacto knife to scrape the debris off.

To Install the new Seals and O-Rings

- 1** Unpack the Tip-Seals. Two tip seals are provided in the kit; one fits the orbiting scroll and the other fits the out-board housing scroll.
- 2** Install the correct tip seal into the groove of the orbiting scroll by:
 - a** Facing the white plastic surface upwards and inserting the foam backing into the groove. Start from the center and work outward.
 - b** Cutting the tip seal off about 1/8" (3 mm) from the outer end of the groove.
- 3** Install the correct tip seal into the groove of the outboard housing scroll as in step 3.
- 4** Place the new main O-ring into the groove in the frame.
- 5** Ensure that the groove is clean.
- 6** Ensure the sealing face of the outboard housing is clean. Carefully replace the outboard housing by lining up the locating pins. Ensure the tip seal has not fallen out of its groove.
- 7** Reinstall (4) M5 bolts (B) and torque the (4) M5 bolts to 5.6 Nm (50 in-lb).
- 8** Unlock the upper part of TPS-compact by sliding the small lever you can find on the bottom side of TPS-compact.
- 9** Reconnect the TPS-compact to the electrical power mains.

Technical Information

- 10 Disable the "TIP-SEAL SERV" warning message by serial interface. To disable the Tip-seal Service message you have to reset the Tip-seal hour-meter (Win 358). Write a "0" on the window 358. The message will be showed again after 5000 h of operation. To check the Tip-seal life you have just to read the Win 358.

Tip Seal Life Reset using A-plus

The screenshot shows a sequence of three windows in a software interface:

- Top window:** A dark blue header with the text "Add new parameter".
- Middle window:** A dark blue header with the text "Enter unknown window". Below the header are three input fields: "Window number:" with the value "538", "Description:" with the value "reset", and "Type:" with radio buttons for "boolean", "numeric" (which is selected), and "string". At the bottom are "OK" and "CANCEL" buttons.
- Bottom window:** A dark blue header with a dropdown arrow and the text "Unknown Windows". Below the header is a list item "reset" with a list icon and an eye icon to its right. Below the list is a highlighted row with the text "538 (reset): reset" and "Value: TIMEOUT". At the bottom is an input field containing "0" and a "WRITE VALUE" button.

To test the Pumping System

- 1 Run the TPS-compact for about 5 seconds. Verify that the front fan is running. If you hear loud noises from the Scroll pump this indicates that the tip seal or main O-ring are possibly out of place.
- 2 Disassemble and repair as necessary.
- 3 The pump is now ready to return to service.

NOTE

Newly installed tip seals may require several hours of run time to seat properly and enable the pump to meet speed and base pressure specifications.

Troubleshooting Chart

Table 8

Problem	Possible Causes	Corrective Actions
TMP doesn't reach NORMAL OPERATION	Scroll Pump motor thermal protector open.	Allow motor to cool. Identify cause of overload.
	Vacuum system leak.	Locate and repair leak.
System will not start	Excessive voltage drop.	Check size and length of cable.
	Fuses burnt-out.	Check fuses and if necessary replace.
Poor ultimate pressure	Vacuum system leak.	Locate and repair leak.
	Water in TPS-compact.	Flush the system with air (use provided IDP3GBKIT) or dry nitrogen.
	Water vapour in TPS-compact.	Flush the system with air (use provided IDP3GBKIT) or dry nitrogen. Install trap or filter.
	Scroll Pump seals worn out.	Replace tip seals.
Hammering noise	Poor conductance to TPS-compact.	Replumb with shorter and/or larger diameter tubing.
	Scroll pump overheated.	Check ventilation to TPS-compact. Check ambient temperature.
Internal Scroll Pump runs intermittently	Debris in Scroll pump.	Flush the TPS-compact. Disassemble pump and inspect.
	Motor thermal protector is cycling open and closed.*	Allow the scroll motor to cool. Identify the cause of the overload.
Overtemperature	Fans clogged or not spinning, problems with turbomolecular pump.	Check if the rotor freely rotates.

- * The internal Scroll Pump is equipped with an auto-reset thermal motor protector. This protector automatically shuts down the pump when it detects an overload condition, and automatically restarts the pump when the motor has cooled to within an acceptable temperature range.

CAUTION

Overtemperature:
To prevent damage to the pump when the temperature exceeds 60 °C, the sensor automatically cuts off the power supply.

Accessories and Spare Parts

Table 9 Accessories and spare parts

TPS-compact	74 FS KF40	74 FS CFF 2.75"	74 FS ISO 63	74 FS CFF 4.5"
115 V	X3583-64000	X3583-64001	X3583-64002	X3583-64003
230 V	X3583-64004	X3583-64005	X3583-64006	X3583-64007
100 V	X3583-64008	X3583-64009	X3583-64010	X3583-64011
Accessories				
DIY Full Range Gauge Kit	9699190	9699199	9699192	9699193
Wheels accessory kit	9698233	9698233	9698233	9698233
Mains cable NEMA plug [3m]	9699958	9699958	9699958	9699958
Mains cable EU plug [3m]	9699957	9699957	9699957	9699957
Inlet screen	9699309	9699328	9699300	9699300
FRG-700 3 m cable	9699960	9699960	9699960	9699960
FRG-700 5 m cable	9699961	9699961	9699961	9699961
Serial cable	9699883	9699883	9699883	9699883
Tip-seal replacement kit	IDP3TS	IDP3TS	IDP3TS	IDP3TS
Vent Valve N.O. 0.5 mm orifice	9699844M005	9699844M005	9699844M005	9699844M005
Mains Cable China Plug	8121-0723	8121-0723	8121-0723	8121-0723
Power cable IEC320 (3 m long)	9499399	9499399	9499399	9499399
Breather vent with adapter	IDP3GBKIT	IDP3GBKIT	IDP3GBKIT	IDP3GBKIT

Technical Information

Table 9 Accessories and spare parts (continued)

TPS-compact	305 FS ISO 100	305 FS CFF 6"	305 FS ISO 160	305 FS CFF 8"
115 V	X3582-64000	X3582-64001	X3582-64002	X3582-64003
230 V	X3582-64004	X3582-64005	X3582-64006	X3582-64007
100 V	X3582-64008	X3582-64009	X3582-64010	X3582-64011
Accessories				
DIY Full Range Gauge Kit	9699194	9699201	9699195	9699202
Wheels accessory kit	9698233	9698233	9698233	9698233
Mains cable NEMA plug [3m]	9699958	9699958	9699958	9699958
Mains cable EU plug [3m]	9699957	9699957	9699957	9699957
Inlet screen	9699302	9699302	9699304	9699304
FRG-700 3 m cable	9699960	9699960	9699960	9699960
FRG-700 5 m cable	9699961	9699961	9699961	9699961
Serial cable	9699883	9699883	9699883	9699883
Tip-seal replacement kit	IDP3TS	IDP3TS	IDP3TS	IDP3TS
Breather vent with adapter	IDP3GBKIT	IDP3GBKIT	IDP3GBKIT	IDP3GBKIT
Vent Valve N.O. 0.5 mm orifice	9699844M005	9699844M005	9699844M005	9699844M005
Mains Cable				
China Plug	8121-0723	8121-0723	8121-0723	8121-0723

For a complete overview of Agilent's extensive product lines, please refer to the Agilent catalog.

Vacuum Products Division

Dear Customer,

Thank you for purchasing an Agilent vacuum product. At Agilent Vacuum Products Division we make every effort to ensure that you will be satisfied with the product and/or service you have purchased.

As part of our Continuous Improvement effort, we ask that you report to us any problem you may have had with the purchase or operation of our products. On the back side you find a Corrective Action request form that you may fill out in the first part and return to us.

This form is intended to supplement normal lines of communications and to resolve problems that existing systems are not addressing in an adequate or timely manner.

Upon receipt of your Corrective Action Request we will determine the Root Cause of the problem and take the necessary actions to eliminate it. You will be contacted by one of our employees who will review the problem with you and update you, with the second part of the same form, on our actions.

Your business is very important to us. Please, take the time and let us know how we can improve.

Sincerely,

A handwritten signature in black ink, appearing to read "Giampaolo LEVI".

Giampaolo LEVI

*Vice President and General Manager
Agilent Vacuum Products Division*

Note: Fax or mail the Customer Request for Action (see backside page) to Agilent Vacuum Products Division (Torino) – Quality Assurance or to your nearest Agilent representative for onward transmission to the same address.

CUSTOMER REQUEST FOR CORRECTIVE / PREVENTIVE / IMPROVEMENT ACTION

TO: AGILENT VACUUM PRODUCTS DIVISION TORINO – QUALITY ASSURANCE FAX

N°: XXXX-011-9979350

ADDRESS: AGILENT TECHNOLOGIES ITALIA S.p.A. – Vacuum Products Division –

Via F.Ili Varian, 54 – 10040 Leini (TO) – Italy

E-MAIL: vpd-qualityassurance_pdl-ext@agilent.com

NAME	COMPANY	FUNCTION
<p>ADDRESS:</p> <p>TEL. N° : FAX N° :</p> <p>E-MAIL:</p>		
<p>PROBLEM / SUGGESTION :</p> <p>REFERENCE INFORMATION (model n°, serial n°, ordering information, time to failure after installation, etc.):</p> <p style="text-align: right;">DATE</p>		
<p>CORRECTIVE ACTION PLAN / ACTUATION (by AGILENT VPD)</p>		<p>LOG N°</p>

XXX = Code for dialing Italy from your country (es. 01139 from USA; 00139 from Japan, etc.)

**Vacuum Products Division
Instructions for returning products**

Dear Customer,

Please follow these instructions whenever one of our products needs to be returned.

Complete the attached **Request for Return form** and send it to Agilent Technologies (see below), taking particular care to include the completed **Health and Safety** declaration Section. No work can be started on your unit until we receive a completed copy of this form.

After evaluating the information, Agilent Technologies will provide you with a **Return Authorization (RA) number** via email or fax, as requested. Note: Depending on the type of return, a Purchase Order may be required at the time the **Request for Return** is submitted. We will quote any necessary services (evaluation, repair, special cleaning, eg).

Product preparation

- Remove all accessories from the core product (e.g. inlet screens, vent valves).
- Prior to shipment and if applicable for your product, drain any oils or other liquids, purge or flush all gasses, and wipe off any excess residue.
- If ordering an Advance Exchange product, please use the packaging from the Advance Exchange to return the defective product.
- Seal the product in a plastic bag, and package product carefully to avoid damage in transit. You are responsible for loss or damage in transit.
- Include a copy of the Health and Safety Declaration in the shipping documentation on the outside of the shipping box of your returning product.
- Clearly label package with RA number. Using the shipping label provided will ensure the proper address and RA number are on the package. Packages shipped to Agilent without a RA clearly written on the outside cannot be accepted and will be returned.
- Return only products for which the RA was issued.

Shipping

- Ship to the location specified on the printable label, which will be sent, along with the RA number, as soon as we have received all of the required information. Customer is responsible for freight charges on returning product.
- Return shipments must comply with all applicable Shipping Regulations (IATA, DOT, ADR, etc.) and carrier requirements.

RETURN THE COMPLETED REQUEST FOR RETURN FORM TO YOUR NEAREST LOCATION:

EUROPE:

Fax: 00 39 011 9979 330
Fax Free: 00 800 345 345 00
Toll Free: 00 800 234 234 00

vpt-customer@agilent.com

NORTH AMERICA:

Fax: 1 781 860 9252
Toll Free: 800 882 7426

vpl-ra@agilent.com

PACIFIC RIM:

please visit our website for individual
office information

<http://www.agilent.com>

TERMS AND CONDITIONS

Please read the terms and conditions below as they apply to all returns and are in addition to the Agilent Technologies Vacuum Product Division – Products and Services Terms of Sale.

- Unless otherwise pre-negotiated, customer is responsible for the freight charges for the returning product. Return shipments must comply with all applicable **Shipping Regulations** (IATA, DOT, etc.) and carrier requirements.
- Agilent Technologies is not responsible for returning customer provided packaging or containers.
- Customers receiving an Advance Exchange product agree to return the defective, rebuildable part to Agilent Technologies **within 15 business days**. Failure to do so, or returning a non-rebuildable part (crashed), will result in an invoice for the non-returned/non-rebuildable part.
- Returns for credit toward the purchase of new or refurbished Products are subject to prior Agilent approval and may incur a restocking fee. Please reference the original purchase order number.
- Units returned for evaluation will be evaluated, and a quote for repair will be issued. If you choose to have the unit repaired, the cost of the evaluation will be deducted from the final repair pricing. A Purchase Order for the final repair price should be issued within 3 weeks of quotation date. Units without a Purchase Order for repair will be returned to the customer, and the evaluation fee will be invoiced.
- Products returned that have not been drained from oil will be disposed.
- A Special Cleaning fee will apply to all exposed products
- If requesting a calibration service, units must be functionally capable of being calibrated.

Vacuum Products Division
Request for Return Form

Customer information
Company :
Address:
Contact Name:
Tel:
Fax:
Email:

Equipment
Product description
Agilent PartNo
Agilent Serial No
Original Purchasing Reference
Failure description
Type of process (for which the equipment was used)

Type of return
Non Billable
Billable
New PO # (hard copy must be submitted with this form):
Exchange
Repair
Upgrade
Consignment/Demo
Calibration
Evaluation
Return for Credit

Health and safety
Substances (please refer to MSDS forms)
The product has been exposed to the following substances:
Trade name
Chemical name
Chemical Symbol
CAS Number
Toxic
Harmful
Corrosive
Reactive
Flammable
Explosive (*)
Radioactive (*)
Biological (*)
Oxidizing
Sensitizer
Other dangerous substances

Goods preparation
If you have replied YES to one of the above questions. Has the product been purged?
If yes, which cleaning agent/method:
Has the product been drained from oil?
I confirm to place this declaration on the outside of the shipping box.

I declare that the above information is true and complete to the best of my knowledge and belief.
I understand and agree to the terms and conditions on page 2 of this document.
Name:
Position:
Date:
Authorized Signature:

NOTE: If a product is received at Agilent which is contaminated with a toxic or hazardous material that was not disclosed, the customer will be held responsible for all costs incurred to ensure the safe handling of the product, and is liable for any harm or injury to Agilent employees as well as to any third party occurring as a result of exposure to toxic or hazardous materials present in the product.

Agilent Vacuum Products Division/Sales and Service Offices

United States

Agilent Technologies

121 Hartwell Avenue
Lexington, MA 02421 - USA
Ph. +1 781 861 7200
Fax: +1 781 860 5437
Toll-Free: +1 800 882 7426
vpl-customer@agilent.com

Netherlands

Agilent Technologies Netherlands B.V.

Customer Contact Center
Laan van Langerhuize 1, toren A-8
1186 DS Amstelveen
Tel. +31 020 547 2600
Fax +31 020 654 5748
customer@agilent.com

Belgium

Agilent Technologies Belgium S.A./N.V.

Customer Contact Center
Pegasus Park
De Kleetlaan 12A bus 12
B-1831 Diegem
Tel. +32 2 404 92 22
Fax +32 2 626 46 30
customer@agilent.com

Brazil

Agilent Technologies Brasil

Avenida Marcos Penteado de Ulhoa
Rodrigues, 939 - 6° andar
Castelo Branco Office Park
Torre Jacarandá - Tamboré
Barueri, Sao Paulo CEP: 06460-040
Toll free: 0800 728 1405

China

Agilent Technologies (China) Co. Ltd

No.3, Wang Jing Bei Lu, Chao Yang District
Beijing, 100102, China
Tel: +86 (0)10 64397888
Fax: +86 (0)10 64392765
Toll free: 400 8206778 (mobile)
Toll free: 800 8206778 (landline)
vacuum.cnmarketing@agilent.com
vpc-customerservice@agilent.com

France

Agilent Technologies
Parc Technopolis - Z.A. de Courtaboeuf
3, avenue du Canada - CS 90263
91978 Les Ulis cedex, France
Tel: +33 (0) 1 64 53 61 15
Fax: +33 (0) 1 64 53 50 01
vpf.sales@agilent.com

Southeast Asia

Agilent Technologies Sales Sdn Bhd

Unit 201, Level 2 uptown 2,
2 Jalan SS21/37, Damansara Uptown
47400 Petaling Jaya, Selangor, Malaysia
Ph. +603 7712 6181
Fax: +603 7727 1239
Toll free: 1 800 880 805
vps-customerservice@agilent.com

India (Sales)

Agilent Technologies India Pvt. Ltd.

Unit Nos 110- 116, & Part of 101 & 109
First Floor, Splendor Forum,
Plot No.-3, District Centre, Jasola
New Delhi-110025
Ph. +91 11 4623 7100
Fax: +91 4623 7105
Toll Free: 1 800 180 1517

Italy

Agilent Technologies Italia S.p.A.

Via F.lli Varian, 54
10040 Leini, (Torino) - Italy
Tel: +39 011 9979 111
Fax: +39 011 9979 350
Toll free: 00 800 234 234 00
vpt-customerservice@agilent.com

Japan

Agilent Technologies Japan, Ltd.

9-1 Takakura-cho Hachioji-city,
Tokyo, Japan
Tel.: +81- 3-5232-1253
Fax: +81-120-565-154
Toll-Free: +81-120-477-111
jp-vvt-sales.pdl-ext@agilent.com

Singapore

Agilent Technologies Singapore Pte. Ltd.

1 Yishun Avenue 7,
Singapore 768923
Tel : (65) 6215 8045
Fax : (65) 6754 0574
Toll free: 1 800 2762622
vps-customerservice@agilent.com

Korea

Agilent Technologies Korea, Ltd.

Ilshin Building 4F
Yongsan-gu Hannam-daero
Seoul Korea 04418
Tel: +82 (0)2 2194 9449
Fax: +82 (0)2 2194 9853
Toll free: 080 222 2452
vpc-customerservice@agilent.com

UK and Ireland

Agilent Technologies LDA UK Limited

Lakeside Cheadle Royal Business Park
Cheadle, Cheshire SK8 3GR,
United Kingdom
Ph. +44 01865291570
Fax +44 01865291571
Toll free: 00 800 234 234 00
Toll free fax: 00 800 345 345 00
vpt-customer@agilent.com

© Agilent Technologies, Inc. 2021
02/2021

India (Service)

Agilent Technologies India Pvt. Ltd.
C-Block, RMZ Centennial Plot Number- 8A, 8B, 8C,
8D,
Doddanakundi Industrial Area, ITPL Road,
Mahadevapura Post, Bangalore- 560048
Ph. +91 80 40614000
Fax: +91 80 40148991

Taiwan

Agilent Technologies Taiwan Limited

No. 20 Gao-shuang Road, Ping-zhen Dist
Tao-Yuan City
32450 Taiwan, R.O.C.
Tel: +886 3 4959004
Toll free: 0800 018 768
vpw-customerservice@agilent.com

Germany and Austria

Agilent Technologies Sales & Services GmbH & Co. KG

Lyoner Str. 20
60 528 Frankfurt am Main
GERMANY
Tel: +49 69 6773 43 2230
Fax: +49 69 6773 43 2250

Mexico

Agilent Technologies

Concepcion Beistegui No 109
Col Del Valle
C.P. 03100 – Mexico, D.F.
Tel.: +52 5 523 9465
Fax: +52 5 523 9472

Other Countries

Agilent Technologies Italia S.p.A.

Via F.lli Varian, 54
10040 Leini, (Torino) - Italy
Tel.: +39 011 997 9111
Fax: +39 011 997 9350
Toll-Free: 00 800 234 234 00
vpt-customer@agilent.com

Customer Support & Service

NORTH AMERICA:

Toll Free: 800 882 7426
vpl-ra@agilent.com
Lexington-service@agilent.com

EUROPE:

Toll Free: 00 800 234 234 00
vpt-customer@agilent.com

PACIFIC RIM:

please visit our website for individual office
information
<http://www.agilent.com/chem/vacuum>

Worldwide Web Site, Catalog and Order On- line:

www.agilent.com/chem/vacuum
Representatives in most countries

In This Book

The manual describes the following:

- Istruzioni per l'uso
- Bedienungshandbuch
- Notice de mode d'emploi
- Manual de instrucciones
- 用户手册
- ユーザーマニュアル
- Instruction for Use
- Technical information

www.agilent.com

© Agilent Technologies, Inc. 2022

Edition L.00, 04/22

87-901-010-01