

Individual Impression Material Units Comparison Chart

INDIVIDUAL IMPRESSION UNITS							
Material	Item Number	Color	Package Contents	Manufacturer	Viscosity	Setting Time	Shore-A Value
Flextime Singles	JKR96C	Pink / Purple	Each	JKR Labs	Low	4 min	30
Gold Velvet II Singles	GV-1	Yellow	Each	JKR Labs	Low	6 min	40
Sihouette Singles	SIL1	Blue	Box of 30	Siemens	Low	3m +/- 20s	26 +/- 2
Dur-A-Sil Equal Singles	RB1	Blue	Each	ipi	Medium	5-6 min	25
Silicone Singles	304	Pink	Box of 100 or by Each	Westone	High	3-5 min	24-30