

New upgrade of electronic equipment

P
e
t
s

Electronic Thermometer

Care for pet health Experience pet love



Medical grade

Fast temperature measurement
Cats and dogs can use it

Oral, axillary and anal can be measured

Three-site measurement

Design of three-position measurement

Anus、Underarm、Oral Cavity Three-site measurement

Satisfy the measuring habits of the owner and pet




Sensitive probe

It takes eight seconds to measure the temperature

Conventional thermometers take 8-10 minutes to reach thermal equilibrium

Pet Can not maintain a posture, easy to cause measurement error




Beep prompt

Prompt for the end of measurement

About 8-10 seconds, the measurement is completed when BB sound is heard

Make it easy for you to master the measurement time




memory function

Record the temperature change of the pet

Digital display, memory function

It is convenient to compare the two measurement data



No hidden danger

Mercury free

No mercury, effectively avoid the hidden danger caused by broken glass or mercury leakage, safe and environmentally friendly

- Mercury is toxic
- Glass is fragile
- Dogs don't like to clip
- Measurement is slow

2000年起, 美国各州开始逐步禁售水银体温计;
2005年, 欧盟开始全面禁止出售水银体温计;
中国政府将在2020年全面禁售水银体温计, 禁用已成为趋势!

Popularization of cat and dog temperature knowledge

The dog's temperature is unstable, Lower in the morning, higher in the afternoon, Just after exercise, body temperature will rise significantly, So we need to record the temperature of the dog during the healthy and quiet period



Dog Body temperature

Cat Body temperature

Small dog: 38.2~38.8°


Cat temperature: 38.0~39.5°

Medium dog: 38.0~38.5°

Big Dog: 37.5~38.5°

ACCURACY UP TO 99%

Accurate / Fast / Convenient / Safe



PET ELECTRONIC THERMOMETER
宠物电子体温计