

# ILLUMINATING LED EYE, FACE & NECK DEVICE & INFUSION TREATMENT

Apply any ampoule of your choice or combine a cocktail to address your skin issues and apply to the skin  
Ampoule choices: Hexapeptide Anti-aging Ampoule, Melanin Ampoule, Mandelic Acid Ampoule, Blue Copper Ampoule, EGF/KGF Ampoule, Filorga NCTF 135 HA, Cytocare 532 Ampoule

## Glycosaminoglycans Conductor Gel or Stop Aging Peptide Essence

Glycosaminoglycans Ingredients

Glycerin, Aloe Barbadensis Leaf Juice, Sorbitol Helianthus Annuus (Sunflower) Seed Oil, Olea Europea (Olive) Oil, Glycosaminoglycans, Sodium Hyaluronate, Panthenol

Glycosaminoglycans (GAGs)—also known as mucopolysaccharides —play key roles in maintaining your health. This has made it a popular ingredient in anti-aging products. The water-binding molecules support skin health, maintaining collagen and elastin as well as hydration.

## ILLUMINATING LED EYE, FACE & NECK DEVICE

Promotes elasticity of the skin

Visible reduction of fine line

Reduces wrinkles

Reduces eye bags

Increases vitality of skin cells

Banish under-eye circles

Eliminates puffiness of Eyes

Better absorption of products by opening pores

Helps with product penetration

Our innovative LED Light Therapy Wand, the ultimate solution for skin rejuvenation, offers a portable design and simple operation. The wand boasts 7 LED colors, each tailored to specific skin needs, from boosting blood circulation to calming inflammation. The wand delivers a heat option for product penetration.

Enjoy a relaxing skincare routine at home, with the 90° rotating massage head catering to every inch of your skin.

Made with high-quality materials and advanced technology, this wand ensures a safe, comfortable, and eco-friendly experience. Don't miss out on this powerful yet gentle skincare tool.

What does the Joli Light Therapy Wand do?

Offers Seven LED Light Therapy:

Red Light Therapy for the face is suitable for all skin types creating blood circulation and collagen production, wrinkles and fine lines.

Green Light Therapy improves pigmentation and melasma

Blue Light Therapy calms and tightens skin

Yellow Light Therapy smooths skin and reduces redness

Purple Light relaxes the skin, improves lymphatic drainage and metabolism,

Cyanine Light soothes, helps with allergies

White Light accelerates tissue metabolism

Galvanic Heat Function:

Heat function is equipped with a 42°C hot compress for soothing warmth, promoting circulation and relieving tension around the eyes and deeper penetration of moisture.

Massages delicate eye area, face and neck and helps to reduce puffiness and dark circles.

Face Lifting helps to lift and tighten the skin around the eyes, promoting a more youthful and firm appearance.

EMS - Microcurrent

Due to the subtle nature of the EMS microcurrent, it is possible that you may not feel a strong sensation while using the product. This is normal and does not indicate a malfunction. We understand that individual experiences may differ, and we appreciate your understanding in this matter.

EMS Electrical Microcurrent Technology utilizes low-level electrical currents that imitate our body's natural bioelectricity. These gentle currents target individual cells, stimulating collagen production and enhancing skin circulation. Consequently, many individuals turn to microcurrent therapy for its anti-aging effects, which include reducing wrinkles and fine lines while promoting better overall skin health.

Galvanic energy is named after Luigi Galvani, an Italian physicist who was studying electric charges in animals in the 18th century, and it's the same type of current that cell phones use. We thank you for that, Mr. Galvani. Its benefits to our skin are numerous: "Galvanic current uses positively and negatively charged ions to push skincare ingredients deeper into the skin, helping to enhance their effectiveness,"

FCC CERTIFICATION ENSURES THAT ELECTRONIC DEVICES COMPLY WITH THE UNITED STATES' STANDARDS FOR RADIOFREQUENCY (RF) AND ELECTROMAGNETIC INTERFERENCE (EMI)<sup>1</sup>. WHEN A PRODUCT IS DEFINED AS "FCC CERTIFIED," IT MEANS IT'S

CE  
VALID SINCE 03/09/2021  
CONFORMITÉ EUROPÉENNE (EUROPEAN CONFORMITY): INDICATES CONFORMITY WITH HEALTH, SAFETY AND ENVIRONMENTAL PROTECTION STANDARDS FOR PRODUCTS SOLD WITHIN THE EUROPEAN ECONOMIC AREA.